



# DiodeVET

*Versatile Veterinary laser for a truly bloodless field*



With 15 years of service, the Diode Vet has established itself as a trustworthy and dependable tool in the clinic.

Designed for minimally invasive, open, or laparoscopic surgery.

Less pain, blood loss, and swelling.

## Why Diode Vet?



### EASY TO USE

*Friendly and intuitive user interface for quick set up. Easily optimize treatment parameters.*



### RELIABLE

*No maintenance or service requirements. No expensive field service calls.*



### COMPACT

*Easily integrate into Doctor's office or operating suite. No hassle of articulated arm.*



**Official training partner of the AIMLA**

**RACE approved**



(845) 277-1700



sales@endoscopy



www.endoscopy.com



## DiodeVET, the trusted name in Veterinary Laser Surgery for over 15 years

*"The diode laser is the most cost effective and efficient wavelength for endoscopic surgery on horses. Contact incision–excision is the forte of the 25W unit. Applicable procedures are epiglottic entrapment, ventriculocordectomy, soft palate scarification, and others."*

– Kenneth E. Sullins, DVM, MS, Diplomate ACVS, Professor of Surgery  
Marion duPont Scott Equine Medical Center

### Specifications:

**Wavelength:** 980nm diode

**Output Power:** 30W Max output

**Dimensions:** 16" x 8.5" x 12"

**Weight:** 18 lbs.

**Handpiece:** Stylus-like

**Operating Modes:** CW or Pulsed (5Hz, 10Hz, or 15 Hz)



**Clinical Applications**