



EVIS EXERA COMPATIBLE GASTROINTESTINAL VIDEOSCOPE

VET TYPE 160-150


IDEAL FOR ROUTINE EXAMINATIONS

With its narrow 8.6 mm diameter insertion tube and wide 2.8 mm diameter channel, this scope offers high image quality, expanded field of view, and great depth of field you need for routine examinations





Main Features for Excellent Performance

- 8.6 mm outer diameter insertion tube and 2.8 mm channel diameter combined with excellent image quality to make this scope ideal for routine examinations.
- Extra-wide 140° field of view enables wider area of site to be easily observed.
- 4-way angulation (210° up, 90° down and 100° right/left) enables complete and comprehensive examination of the upper digestive tract.
- Ergonomically designed grip enhances scope maneuverability while redesigned easy-to-access knobs and user-selectable switches improve operability.
- Superior image quality is delivered even in the expanded full height mode when combined with the dedicated CV-160 video system center.
- Fully compatible with the CV-190/180/160/140/100. Scope ID function stores individual scope information and displays it on the monitor to facilitate endoscopy suite management.



Specifications

Optical System	Field of view	140°
	Direction of view	0° Forward viewing
	Depth of field	3 to 100 mm
Distal End	Outer diameter	8.6 mm
Insertion Tube	Outer diameter	8.6 mm
Bending Section	Angulation range	Up 210°, Down 90°, Right 100°, Left 100°
Working Length		1500 mm
Total Length		1815 mm
Instrument Channel	Inner diameter	2.8 mm
	Minimum visible distance	3 mm from the distal end
	Endo-Therapy accessory entrance/exit position in field of view	



Endoscopy Support Services, Inc.
 3 Fallsview Lane
 Brewster, NY 10509
 845-277-1700
 845-277-7300 Fax
www.endoscopy.com