

Metricide™ High-Level Disinfectant

Metricide OPA Plus is a high level disinfectant for reprocessing heat sensitive semi-critical medical devices for which sterilization is not suitable.

1 gallon container: \$60.00
Case of 4x 1 gal. containers: \$240.00
Test strips (2x bottles of 100): \$195.00



PUREChannel™ Flush

PureChannel™ Flush Instrument Channel Cleaner is a powerful surfactant-enzyme formulation designed for soaking and flushing channeled devices.
Sold as one (1) 32oz bottle.



Part/SKU#: PCF-32OZ : \$35.00 ea.



SKU#: NC-490908

Veterinary COMPANION \$21⁹⁵
Disinfectant Wipes (160 ct.) Per Container



- Effectively kills Parvovirus and Feline Calicivirus
- Provides effective cleaning and disinfection in a single operation
- Fresh citrus scent
- Non-corrosive even after repetitive uses
- No rinsing required
- Improve patient turnover
- Staff-friendly

Disposable Channel Cleaning Brushes

DCB-1.72306 bi-directional for 2.0-4.5mm channels (6mm brush, 1.7mm x 230cm) Pkg of 10.....\$45.00

Best Sellers - Best Value!

DCB-1.01202 bi-directional for 1.2-1.8mm channels (2mm brush, 1.0mm x 120cm) Pkg of 10.....\$45.00

DTP-02 Scope Sponge

Distal Tip Protection & Cleaning Sponge



Shaped sponge designed to clean the insertion tubes of flexible endoscopes. Can also be used as distal tip protection padding during transport.


*Each Sponge: \$3.00

*Also sold as a fifty (50) pack for: \$75.00

Don't have the cash for new?

We have a great selection of used/refurbished equipment to fit you needs and budget.

endoscopy.com/used



ESS ATT-Series Airtight Leak Testers



OLYMPUS • PENTAX • WOLF • OPTIM • STORZ • FUJINON

- Hand-held design
- Easy to use
- Helps you protect your endoscopes

Price each: \$225⁰⁰

Universal Camera Drapes



- Drapes come sterilized
- Reduces turnaround time between cases
- Eliminates camera breakdowns caused by repeated soakings
- Reduces repair costs due to sterilant damage
- Eliminates cost of sterilants and exposure to harmful vapors
- Specifications: 7 in x 96 in.

Box of 50: \$395⁰⁰

(Part/SKU: ESS-6242 CAMDRAPES)

Find more cleaning products on our website: endoscopy.com/clean

Eyepiece Adapters

Adapters to allow rigid camera endocouplers and the ClearSCOPE™ Smartphone Adaptor to connect to your flexible fiberoptic endoscope eyepiece.



Olympus® OES Adapter
Fits on Olympus OES eyepieces. Use set screw to hold adapter.

Part/SKU#: WA-01
Price: \$225.00 ea.



Olympus® pre-OES Adapter
Fits on older Olympus pre-OES eyepieces.

Part/SKU#: WA-02
Discount Price: \$50.00 ea.



Pentax® Adapter
Fits on Pentax eyepieces. Use set screw to hold adapter.

Part/SKU#: WA-03
Price: \$225.00 ea.

OLYMPUS®

endoscopy.com/OlympusVet

ENDOSCOPES

COMPLETE SYSTEMS



PROCESSORS/
LIGHT SOURCES



- Certified Pre-Owned - Direct from the Manufacturer
- High-Resolution HDTV
- Narrow-Band Imaging (NBI)
- Electronic Magnification
- Olympus Quality & Reliability
- Included Olympus-brand accessories
- One-year Olympus Manufacturer's Warranty

ADD RIGID ENDOSCOPES TO YOUR SYSTEM

ESS Veterinary Otoscope Set



- Zero (0) deg. Otoscope w/2.0mm working channel
- Accepts 5FR (1.67mm) instruments
- Dual stopcock valve with luer connector, two (2) channel caps
- Disposable cleaning brush
- Padded storage case

EndoTech™ Veterinary Video Endoscopy System



- Crystal-clear image quality
- Bright, long-lasting LED lighting
- Lightweight and portable
- Easy to use
- Video and image still capture

endoscopy.com/EndoTech

TEAC UR-4MD HD Video Recorder



Medical Grade HD Video Recorder with auto sensing video connections and automatic video resolution detection. Saves to USB memory/thumb drives. NTSC and PAL compatibility.

endoscopy.com/TEAC



endoscopy support services, inc.

best
Your source for everything endoscopy since 1989™

VETERINARY MINI CATALOG

We're the Exclusive Veterinary Distributor of:

OLYMPUS®

Endoscopy Products



1-800-DVM-ENDO

For more of our products, visit us on the web www.endoscopy.com

VETERINARY ENDOSCOPY INSTRUMENTS - BIOPSIES, RETRIEVAL, REPLACEMENTS AND MORE!

REUSABLE • REPAIRABLE • AUTOCLAVABLE

- Stainless steel construction
- Ergonomically designed handle
- Sterilize via autoclave, ETO gas or soaking methods

Price: \$290.00 each
90 Day Warranty

***Standard short oval cup:**
 QFC-3015: 1.8mm x 35cm
 QFC-1015: 1.8mm x 120cm
 QFC-1016: 1.8mm x 175cm
 QFC-1096: 1.8mm x 190cm
 QFC-1216: 2.4mm x 120cm
 QFC-1266: 2.4mm x 170cm
 QFC-1296: 2.4mm x 225cm

***Long ellipsoid cup:**
 QFC-3014: 1.8mm x 35cm
 QFC-1027: 1.8mm x 120cm
 QFC-1010: 1.8mm x 175cm
 QFC-1092: 1.8mm x 190cm
 QFC-1214: 2.4mm x 120cm
 QFC-1212: 2.4mm x 170cm
 QFC-1292: 2.4mm x 225cm

***Long ellipsoid cup w/needle:**
 QFC-3019: 1.8mm x 35cm
 QFC-1022: 1.8mm x 120cm
 QFC-1012: 1.8mm x 175cm
 QFC-1090: 1.8mm x 190cm
 QFC-1211: 2.4mm x 120cm
 QFC-1210: 2.4mm x 170cm
 QFC-1290: 2.4mm x 225cm

***Serrated cup:**
 QFC-3013: 1.8mm x 35cm
 QFC-1024: 1.8mm x 120cm
 QFC-1064: 1.8mm x 175cm
 QFC-1094: 1.8mm x 190cm
 QFC-1234: 2.4mm x 120cm
 QFC-1264: 2.4mm x 170cm
 QFC-1294: 2.4mm x 225cm

***Serrated cup with needle:**
 QFC-1087: 1.8mm x 120cm
 QFC-1077: 1.8mm x 175cm
 QFC-1097: 1.8mm x 190cm
 QFC-1237: 2.4mm x 120cm
 QFC-1277: 2.4mm x 170cm
 QFC-1297: 2.4mm x 225cm

***Alligator jaw retriever:**
 QGR-3016: 1.8mm x 35cm
 QGR-1023: 1.8mm x 120cm
 QGR-1063: 1.8mm x 175cm
 QGR-1093: 1.8mm x 190cm
 QGR-1233: 2.4mm x 120cm
 QGR-1263: 2.4mm x 170cm
 QGR-1293: 2.4mm x 225cm

***Rat tooth retriever:**
 QGR-3017: 1.8mm x 35cm
 QGR-1025: 1.8mm x 120cm
 QGR-1065: 1.8mm x 175cm
 QGR-1095: 1.8mm x 190cm
 QGR-1235: 2.4mm x 120cm
 QGR-1265: 2.4mm x 170cm
 QGR-1295: 2.4mm x 225cm

Three-prong grasper:
 QGR-3018: 1.8mm x 35cm
 QGR-1029: 1.8mm x 120cm
 QGR-1061: 1.8mm x 175cm
 QGR-1091: 1.8mm x 190cm
 QGR-1231: 2.4mm x 120cm
 QGR-1261: 2.4mm x 170cm
 QGR-1291: 2.4mm x 225cm

Forked jaw grasper:
 QGR-3009: 1.8mm x 35cm
 QGR-1026: 1.8mm x 120cm
 QGR-1060: 1.8mm x 175cm
 QGR-1080: 1.8mm x 190cm
 QGR-1226: 2.4mm x 120cm
 QGR-1260: 2.4mm x 170cm
 QGR-1280: 2.4mm x 225cm

Tripod retriever:
 QGR-1032: 1.8mm x 120cm
 QGR-1062: 1.8mm x 175cm
 QGR-1098: 1.8mm x 190cm
 QGR-1232: 2.4mm x 120cm
 QGR-1262: 2.4mm x 170cm
 QGR-1298: 2.4mm x 225cm

Silicone Lubrication Oils

SIL-01
 3x Tubes for forceps & graspers lubrication
\$15.00/pkg of 3

SIL-03B
 15ml Bottle for endoscope button O-rings & seals
\$38.00 ea.

Olympus®-Compatible Buttons

• HU-AW196 (Oly MB-196 Compatible)
 Air/water valve for 10,20,30,100 & 130 series scopes

• HU-S197 (Oly MB-197 Compatible)
 Suction valve for 10,20,30, OSF-3, OSF-V6, EVIS 100/130 series scopes

• HU-AW438 (Oly MH-438 Compatible)
 Air/water valve for 140 series scopes

• HU-S443 (Oly MH-443 Compatible)
 suction valve for 140 series scopes

All Buttons: \$125.00 ea.

Replacement Water Bottles

WB-140 pictured

WB-PEN - \$250.00
 (non-autoclavable for Pentax)

WB-100 - \$175.00
 (non-autoclavable for Olympus 100-130 series)

WB-140 - \$295.00
 (non-autoclavable for Olympus 140 series)



NET-STYLE RETRIEVERS

F012 (Roth Net): 1.8mm x 160cm.....\$160.00
 PTW-2323-30 (Vet Net): 2.3mm x 230cm.....\$160.00

OVAL LOOP RETRIEVAL SNARES

SD-T-18 25200S: 1.8mm x 200cm.....\$30.00
 SD-T-24 25230S: 2.4mm x 230cm.....\$30.00

HEXAGON LOOP POLYPECTOMY SNARE

SD-L-24 25230S (2.4mm x 230cm; 25mm loop)...\$30.00
 SD-L-24 25230 (2.4mm x 230cm; 25mm loop; Bx/10)..\$200.00
 SD-L-18 25200 (1.8mm x 200cm; 25mm loop; Bx/10)..\$200.00

OVAL LOOP GRABBER SNARES

F026: 1.8mm x 240cm\$40.00
 F027: 2.3mm x 240cm\$40.00

JUMBO OVAL LOOP RETRIEVAL SNARES

F028: 1.8mm x 240cm.....\$40.00
 F029: 2.4mm x 240cm.....\$40.00

JUMBO OVAL LOOP GRABBER SNARES

F048: 1.8mm x 240cm.....\$50.00
 F049: 2.3mm x 240cm.....\$50.00

THREE PRONG GRASPING FORCEPS

F024: 1.8mm x 240cm.....\$135.00
 F025: 2.3mm x 240cm.....\$135.00

FOUR PRONG GRASPING FORCEPS

F014: 1.8mm x 240cm.....\$135.00
 F015: 2.3mm x 240cm.....\$135.00

CYTOLOGY BRUSHES

BC-L-2 18180-S: 1.8mm x 180cm.....\$17.00
 (box of 10: \$150.00)
 BC-L-1 24230-S: 2.4mm x 230cm.....\$17.00
 (box of 10: \$150.00)

ASPIRATION CATHETERS

ASPC0: 1.8mm x 30cm (for otoscopes).....\$17.00
 PW-B-1820-S: 1.8mm x 200cm.....\$17.00
 PW-B-2423-S: 2.4mm x 230cm.....\$17.00
 ASPC3: 2.3mm x 480cm (cut to size).....\$40.00

BIOPSY TRANSFER CASSETTES

BTC-1: Package of 25.....\$35.00

Veterinary Basket Retrievers

Stainless Steel Four Wire Basket Retriever
 • F020: 1.8mm x 240cm
 • F021: 2.3mm x 240cm

Stainless Steel Helix Basket Retriever
 • F040: 1.8mm x 180cm
 • F041: 2.3mm x 240cm

"Steeple™" Stainless Steel Jumbo 4-Wire Basket Retriever
 • F050: 1.8mm x 240cm
 • F051: 2.3mm x 240cm

Nitinol Helix-Style Basket Retriever
 • F046: 1.8mm x 180cm
 • F047: 2.4mm x 230cm

Nitinol Angular Basket Retriever
 • F042: 1.8mm x 180cm
 • F043: 2.4mm x 230cm

Nitinol Ball-Style Basket Retriever
 • F044: 1.8mm x 180cm
 • F045: 2.4mm x 230cm

Steeple™ Price: \$150.00 ea.

Nitinol retrievers: \$130.00 ea.

"The Ratigator"

Combo rat tooth & alligator grasper
AUTOCLAVABLE • REPAIRABLE

- TMS-1821RG - 1.8mm x 210cm
- TMS-2423RG - 2.4mm x 240cm

Price each: \$290.00

"The Pelican"

Budget-Priced Combo rat tooth & alligator grasper
NON-AUTOCLAVABLE • SEMI-DISPOSABLE

- FG-AR-1SL 23230C* (Pelican Grasper, 2.3mm x 230cm)
- FG-AR-1SS 23230C* (Mini Pelican Grasper, 2.3mm x 230cm)

Price each: \$125.00

*Disposable Retrievers

- FG-A-1S 23230C-S Alligator jaw graspers, 2.3mm x 230cm
- FG-P-3 18200-S Tripod graspers, 1.8mm x 200cm
- FG-P-3 24230-S Tripod graspers, 2.4mm x 230cm

\$125.00 ea./10 for \$700.00

**"Disposable": use until broken/unusable, then dispose of instrument

* Custom length of 350cm: \$395.00 ea. (NOTE: special order item, no returns/refunds, warranty-only)

Veterinary Disposable Forceps

Box of 10 for \$175.00
NON-AUTOCLAVABLE • DISPOSABLE

- FB-C-1S 23230C Oval cup; 2.3mm x 230cm
- FB-C-1ST 23230C Serrated cup; 2.3mm x 230cm
- FB-C-2 18200C Oval cup; 1.8mm x 200cm
- FB-D-1S 23230SC Oval cup w/spike; 2.3mm x 230cm
- FB-D-1ST 23230 Serrated cup w/spike; 2.3mm x 230cm

Disposable Biopsy Valves

Part/SKU#: DBV-1
 Biopsy valves for Olympus and Fujinon Gastroscopes, Colonoscopes and other larger-channel endoscopes

Part/SKU#: DBV-2
 Biopsy valves for Pentax flexible endoscopes

Package of 10: \$50.00

Part/SKU#: DBV-3
 Biopsy valves for Olympus

Package of 10: \$250.00

ESS MLS1 LED Light Source
COMPACT • BATTERY-POWERED

Portable, bright LED light source for rigid or flexible endoscopes. Runs on two (2) Li-Ion rechargeable batteries (included).

\$995.00 ea.

MLS1 MLS1F

ESS MLS1 LED Light Source

COMPACT • BATTERY-POWERED

Portable, bright LED light source for rigid or flexible endoscopes. Runs on two (2) Li-Ion rechargeable batteries (included).

\$995.00 ea.

*International plug adapter power pack available. Request when ordering.

Endoscopy Storage Systems™

SS-1 ScopeSecure™
 • Holds one (1) endoscope and two (2) forceps
 • Size: 9.5"W x 4.19"H x 4.125"D
 Price: \$89.00 ea.

SS-2 ScopeSecure™ II
 • Holds two (2) endoscopes
 • Size: 14.0"W x 4.31"H x 1.5"D
 Price: \$89.00 ea.

SS-M ScopeSecure™ modular system
 • Holds one (1) endoscope per rack
 • Purchase multiple racks to design a multi-scope storage system
 • Size: 3.75"W x 4.31"H x 5.125"D
 Price: \$89.00 ea.*
*Purchase one at regular price, each additional: \$45.00 ea.

SS-3 ScopeSecure™ III
 • Holds one (1) ENT endoscope, bronchoscope, nasopharyngoscope, and other small flexible scopes
 • Size: 8.5"W x 4.2"H x 4.2"D
 Price: \$58.00 ea.

FS-1 ForcepSecure™
 • Holds six (6) flexible forceps/graspers
 • Size: 10.0"W x 4.25"H x 5.188"D
 Price: \$69.00 ea.

ISS-1 Instrument storage system
 Wall-mounted peg rack with snap-close storage bags to store your flexible endoscope forceps & graspers
 Price: \$29.95 ea.
 ISS-RB: Set of five (5) replacement bags: \$14.95/set

SPECIFICATIONS:

- Custom wall rack with mounting hardware and instructions (5.75"W x 2.75"H x 4.75"D)
- Six (6) vented snap-close hanging storage bags (12.25"H x 9.25"W)
 - 5 bags for instrument storage
 - 1 set of illustrated instrument ID labels (affix to bag for easy identification of contents)
 - 1 bag for broken instrument(s) to be sent for repair
 - ESS "Speed-Up" repair form & three (3) "Out For Repair" cards
- **Three (3) \$10.00 discount coupons to use for subsequent purchases**

ATT-H Airtight tester holder
 • Designed specifically for the ATT Series testers
 • Made of durable plastic
 • Easily cleaned with soap and water
 • Mountable on a wall or endoscopy cart
 • Mounting hardware and instructions included
 Price: \$24.95 ea.

www.endoscopy.com