



**endoscopy
support
services, inc.**



best
Your source for everything endoscopy since 1989™

3 Fallsview Lane • Brewster, NY 10509
Ph: (845) 277-1700 • Fx: (845) 277-7300

Web: www.endoscopy.com

G.P. Endoscopy Made Easy



**2-Way, Rotating
Veterinary
Bronchoscope/
General Use
Videoscopy**

MODEL No.	DIAMETER	WORKING LENGTH	BIOPSY CHANNEL	BENDING	ROTATION	RESOLUTION
EU-120	2.8mm	610mm	1.2mm	Up: 210° Down: 130°	Left/Right: 120°	160P (SD)
EU-130	3.8mm	1000mm	1.5mm	Up: 180° Down: 130°	Left/Right: 120°	160P (SD)
EU-140H	4.8mm	1000mm	2.0mm	Up: 180° Down: 130°	Left/Right: 120°	720P (HD)
EU-150H	5.8mm	1500mm	2.8mm	Up: 180° Down: 130°	Left/Right: 120°	720P (HD)

MODEL No.	DIAMETER	WORKING LENGTH	BIOPSY CHANNEL	RANGE OF BENDING		VIEW ANGLE	RESOLUTION
GU-160H	6.8mm	150cm	2.2mm	Up: 210° Down: 130°	Left/Right: 100°	120°	720P
GU-170HQ	7.8mm	150cm	2.8mm	Up: 180° Down: 130°	Left/Right: 100°	120°	720P
GU-190FHP	9.2mm	150cm	2.8mm	Up: 180° Down: 130°	Left/Right: 100°	140°	1080P
GU-190FHPL	9.2mm	300cm	2.8mm	Up: 180° Down: 130°	Left/Right: 100°	140°	1080P
GU-1B0H-350	11.5mm	350cm	3.2mm	Up: 180° Down: 130°	Left/Right: 160°	140°	1080P

**4-Way Veterinary
G.I./Colono
General Use
Videoscopy**



add the

**VHP-1315
All-In-One Portable**

- Battery-operated
- 1920 x 1080p image resolution
- Works like a Windows® PC with iCap® v5* pre-installed



add an
**iCap® HD USB
Camera**

1080p High Definition
video quality

**iCap® USB Camera
works with virtually
all Rigid Scopes**

HD
High Definition

*Upgrade from iCap® v5 to DICOM v8 for additional cost. Ask your ESS Sales Consultant for details.



endoscopy.com/VetUSB



USB Disposable Video Bronchoscope



- High-Definition video
- 2.8mm insertion tube diameter
- 60cm (600mm) working length
- 1.2mm biopsy channel
- Bending: 230° Up, 180° down; 120° L/R Rotation
- Model number: BRD-1PF
- Comes as two (2) Sterile-Sealed, ready-to-use videoscopes with USB adapter, biopsy valve and suction valve



endoscopy.com/VetUSB

HD Disposable Video Bronchoscopy System



15.6" HD Touchscreen Veterinary Video Processor Model#: MS-8



- Compact, portable easy-to-use design
- Videoscope with biopsy channel & 2-Way deflection (210 degrees Up & Down)
- 5.0mm insertion tube outer diameter
- 100cm (1-meter) working length
- 2.2mm biopsy channel
- 'Chip-On-Tip' CMOS sensor & bright LEDs at the tip
- HD imaging on battery-powered processor at 1024x768p (Attach an external HDMI monitor for 1920x1080p viewing!)
- Comes pre-sterilized - *use right out of the box!*

endoscopy.com/VetBronch



ESS Kit: Ambu® Vet E/R Disposable Foreign Body Retrieval Set (Part/SKU#: ESS-AmbuER1)

Ambu 5.0mm Disposable Endoscope Set with 2x Ambu 5.0mm aScopes, a View Monitor/Processor, TMS-1821 Ratigator, loop snare, and cleaning brushes. Use in veterinary emergency rooms for quick foreign body retrieval in small animals. Leasing available.

endoscopy.com/Ambu



<p>Elevated Comfort Providing the best experience for patients by reducing injection pain.</p>	<p>One-Handed Use The only device to offer single-handed ability to numb and inject quickly.</p>	<p>Fast-Acting Coolject mist has the quickest numbing onset & time to injection only takes 4-10 seconds.</p>
<p>Targeted Application Targeted blanched area for immediate, accurate needle penetration.</p>	<p>Cost-Efficient Coolject provides over twice the amount of uses per ounce than the competition.</p>	<p>Non-Toxic A safer alternative to ethyl chloride, is hazmat free, and non-flammable.</p>



endoscopy.com/VapoCool

SonoScape

Veterinary Video Endoscopy System



SonoScape™ System
(Cart, Monitor, V-2000, V-10M Videoscope & Opt. Camera Head)



V-Series Scopes

- Hi-Definition video image
- Easy, single-connect
- 4-way articulation
- Air/Water & Suction capable
- On-scope operation buttons
- 1-meter x 6.7mm, 2.0mm Chan.
- 1.5M x 11mm; 2.8mm Chan.



V-2000 Processor/Light Source


- Bright LED lighting system
- Touch-screen operation
- Single-connection for Videoscopes
- Optional Rigid Camera attachment
- Video & Still image capture to USB memory

endoscopy.com/SonoVet




SKU#: SSV-V-10S

V-10S Videoscope: 5.4mm Insertion Tube diameter, 110cm (1100mm) Working Length, 2.2mm Biopsy Channel (use 1.8mm diameter instruments)



SKU#: SSV-V-10M

V-10M Videoscope: 9.0mm Insertion Tube diameter, 150cm (1500mm) Working Length, 2.8mm Biopsy Channel (use up to 2.5mm diameter instruments)



SKU#: SSV-V-10L

V-10L Large Animal Videoscope: 12.0mm Insertion Tube diameter, 330cm (3300mm) Working Length, 3.8mm Biopsy Channel (use up to 3.4mm diameter instruments)



SonoScape Camera Head with Focal Lens (19-38mm)

10-110-201V Camera head with focal lens 19-38mm - 60fps full HD image, CMOS image sensor, Standard C-Mount (rigid eyecup endoscopes).

Attaches to the SonoScape V-2000 processor. Requires a separate light guide for light delivery (sold separately).

SKU#:
SSV-10-110-201V

AVC-1 HD Video Endoscopy System



HD Image Sensor High-Quality Optics



Digital Zoom up to 2x



Hemoglobin Enhancement (HbE)



140° Field of View



LCD Touch Screen Operation



AOHUA




endoscopy.com/AOHUA



EVIS EXERA III



NARROW-BAND IMAGING



190 HD Endoscopy System

The 190 Series EVIS EXERA III delivers high-definition imaging for the best view of your subject. **Narrow-Band Imaging (NBI)** highlights blood vessels so you can see the most minute irritations and bleeding areas. The CLV-190 light source is fully connected to the CV-190 processor so the need for a video connecting "pigtail" with 190-series endoscopes is no longer needed. CV-190 had a connection slot for a video connector to use pre-190 Olympus videoscopes with the system.



EVIS EXERA III COLONOVIDEOSCOPE PCF-PH190L/I

The EVIS EXERA III PCF-PH190L colonovideoscope delivers **HDTV** image quality and **Narrow Band Imaging (NBI)** for enhanced visualization in an ultra-slim scope. An outer diameter of only 9.7 mm along with **Responsive Insertion Technology (RIT)** helps ease insertion.




EVIS EXERA III GASTROINTESTINAL VIDEOSCOPE GIF-H190

Olympus GIF-H190 EVIS EXERA III High-Definition, Narrow Band Imaging compatible, video gastroscope with full screen HDTV image, forward water jet function. Specifications: 140 degree field of view; Angulation: Up 210 degree, Down 90 degree, Right and left, 100 degree Outer Diameter: 9.2mm; Working Length: 103cm; Instrument Channel Diameter: 2.8mm



EVIS EXERA III GASTROINTESTINAL VIDEOSCOPE GIF-XP190N

The EVIS EXERA III GIF-XP190N Veterinary gastroscope delivers excellent image quality and illumination, enhanced **NBI** capabilities, and a wide 140° field of view. Its 5.4 mm outer diameter makes the GIF-XP190N a powerful tool for diagnostic endoscopy whenever an ultra-slim scope is needed.



URETERO-RENO VIDEOSCOPE URF-V2

The URF-V2 achieves a new standard in smooth insertability and high image quality. As slim as a fiberscope, yet with the image quality of a videoscope, every feature contributes to this scope's versatile performance. Improved visualization, facilitated maneuverability and reduced operator stress are now within reach.


OTV-S7Pro-H-HD-12E Camera Head

Key Benefits

- **High resolution images:** With 1080i HD output, the camera head provides images for incredibly sharp, accurate imaging.
- **Full autoclavability:** The camera head is tightly sealed, preventing lens fogging even after autoclave steam sterilization reprocessing.
- **Programmable buttons:** The control section has three built-in switches, allowing programming of frequently used functions, such as electronic zoom and image capture for convenient remote control.
- **One touch zoom:** The image can be enlarged to 1.2 or 1.5 times the original size, and high-definition, large images are maintained even in zoom mode.
- **Quick-lock connection:** Q models feature a quick-lock connection, allowing the camera head to be snap locked onto the telescope to ensure a tight fit.



Connects to the VISERA OTV-S7Pro, CV-180 and CV-190 Processors.



URETERO-RENO VIDEOSCOPE URF TYPE V



■ URF-V Specifications

Optical System	Field of View	90°
	Direction of View	0° (Forward)
	Depth of Field	2 ~ 50 mm
Insertion Tube	Distal End Outer Diameter	8.5 Fr.
	Insertion Tube Outer Diameter	9.9 Fr. (3.3 mm)
	Working Length	670 mm
Instrument Channel	Inner Channel Diameter	3.6 Fr. (1.2 mm)
Bending Section	Angulation Range	Up 180°/Down 275°
Total Length		980 mm



Popular Veterinary Instruments & Accessories



Veterinary Disposable Forceps Box of 10 for \$200⁰⁰

NON-AUTOCCLAVABLE - DISPOSABLE

- FB-C-IS 23230C Oval cup; 2.3mm x 230cm
- FB-C-IST 23230C Serrated cup; 2.3mm x 230cm
- FB-C-2 18200C Oval cup; 1.8mm x 200cm
- FB-C-2T 18200C Serrated cup; 1.8mm x 200cm
- FB-D-IS 23230C Oval cup w/spike; 2.3mm x 230cm
- FB-D-IST 23230 Serrated cup w/spike; 2.3mm x 230cm



"Steeple"[™] Stainless Steel Jumbo 4-Wire Basket Retriever

- F050: 1.8mm x 240cm
- F051: 2.3mm x 240cm



Steeple Price:
**\$198⁰⁰
ea.**



"The Ratigator"[™] Combo rat tooth & alligator grasper **AUTOCCLAVABLE - REPAIRABLE**

- TMS-1821RG - 1.8mm x 210cm
- TMS-2423RG - 2.4mm x 240cm

Price each: **\$340.00**



JUMBO OVAL LOOP RETRIEVAL SNARES

- F028: 1.8mm x 240cm - \$65.00
- F029: 2.4mm x 240cm - \$65.00


JUMBO OVAL LOOP GRABBER SNARES

- F048: 1.8mm x 240cm - \$65.00
- F049: 2.3mm x 240cm - \$65.00



CYTOLOGY BRUSHES

- BC-L-2 18180-S: 1.8mm x 180cm - \$27.00
(box of 10: \$243.00)
- BC-L-1 24230-S: 2.4mm x 230cm - \$27.00
(box of 10: \$243.00)



NET-STYLE RETRIEVERS

- F012 (Roth Net): 1.8mm x 160cm.....\$165.00
- PTW-2323-30 (Vet Net): 2.3mm x 230cm.....\$165.00

ESS ATT-Series Airtight Leak Testers



- Hand-held design
- Easy to use
- Helps you protect your endoscopes

Price each:
\$250⁰⁰



ASPIRATION CATHETERS

- ASPC0: 1.8mm x 30cm (for otoscopes).....\$27.00
- PW-B-1820-S: 1.8mm x 200cm.....\$22.00
- PW-B-2423-S: 2.4mm x 230cm.....\$22.00
- ASPC3: 2.3mm x 480cm (cut to size).....\$50.00



Neogen Veterinary
COMPANION Disinfectant
Wipes 160 ct. Pre-moistened
towelettes, 6in x 8in (2lbs,
4oz)



- Effectively kills Parvovirus and Feline Calicivirus
- Provides effective cleaning, disinfection and odor control in a single step
- Kills 99% of germs
- Non-corrosive even after repetitive uses
- No rinsing required
- Helps improve patient turnover efficiency in clinic environment

SKU #: NC-490908
Price ea.: \$21.95

Universal Camera Drapes



- Drapes come sterilized
- Reduces turnaround time between cases
- Eliminates camera breakdowns caused by repeated soakings
- Reduces repair costs due to sterilit damage
- Eliminates cost of sterilants and exposure to harmful vapors
- Specifications: 7 in x 96 in.

Box of 50: **\$395⁰⁰**

(Part/SKU: ESS-6242 CAMDRAPES)



Find more cleaning products on
our website: endoscopy.com/clean

endoscopy.com



SurgASSIST[™]

INSUFFLATOR & ELECTROSURGICAL GENERATOR

Streamline your operating room with the **SurgASSIST[™] CO₂ Insufflator and Electrosurgical Generator Combination**, an innovative solution that brings CO₂ insufflation and electrosurgical capabilities together in a single, compact system.

Designed to enhance efficiency and save valuable space, the SurgASSIST[™] delivers unmatched functionality with Biovision's hallmark ease of use. Its insufflation system is tailored for laparoscopic procedures, featuring crucial tools like the over- and under-pressure level indicator for safety and a Patient Pressure Select switch that allows quick adjustments to operating CO₂ pressures. The integrated electrosurgical generator offers precise power adjustment and reliable performance for cutting and coagulating, all conveniently controlled through foot-switch activation. Whether performing everyday clinical tasks or advanced minimally invasive procedures, the SurgASSIST[™] is engineered to meet your needs with simplicity, reliability, and maximum functionality.

Upgrade your clinical practice with this all-in-one solution and bring efficiency and precision to every procedure.

endoscopy.com/SurgASSIST




SPiTZ[®]

VET CARE

- Oscopes, Cystoscopes, Ureterscopes
- Arthroscopes, Sheaths & Instruments
- Veterinary Shaver System
- Arthroscopy Pumps
- Laparoscopes, Sheaths & Instruments
- Veterinary Insufflators
- HD/4K Video Camera Systems
- LED Daylight-Quality Light Sources

and much more!



endoscopy.com/SPiTZ

Rigid Endoscopy



Vetina
CS3 HD Video
System



ARTHROSCOPY



**OTO/RHINO/
LARYNGOSCOPY**



CYSTOSCOPY



LAPAROSCOPY

endoscopy.com/mindray

MICROLINE
SURGICAL



**MiFusion®
MiSeal
Vessel
Sealing
System**

Disposable handles
with connection cords
available in
3 different lengths.

The MiSeal Thermal Ligating System offers surgical instruments used to seal and divide tissue utilizing proprietary Thermal Fusion technology for Veterinarians.

endoscopy.com/MiSeal

**RICHARD
WOLF**



POWERBLADE
Bipolar Cutting & Coagulation System



PATENTED BIPOLAR
CUTTING BLADE

For Use In
**Laparoscopic
Surgical Procedures**

Contact for pricing & availability: sales@endoscopy.com

CLOUD ENDO

**REVOLUTIONIZING
MEDICAL VIDEO & IMAGING**
The Missing Link to PACS, DICOM, & Medical Video

www.cloudendo.com



Streamline Workflow

CloudEndo™ goes beyond image management with our integrated scheduling module and comprehensive Modality Worklist (MWL) functionality.

Seamlessly schedule appointments and effortlessly retrieve patient data, enabling a more efficient and organized workflow for your practice.

BENEFITS

- Smooth onboarding
- User-friendly interface
- Effortless data access
- High-quality video streaming
- Streamlined reporting
- Secure data storage
- Collaborative sharing
- Reliable support

Send DiCOM Imaging from:


- Video Endoscopy
- MRI & CT
- Digital X-Ray
- Video Microscopy
- Ultrasound

Share Securely to anyone's device:

- Desktop Computer
- Laptop Computer
- Tablet & Smartphone



- Capture still images and video
- **Writes DiCOM information for both!**
- Annotate images & generate reports
- **One-Click CloudEndo™ Archiving**



endoscopy.com/iCap8

Service Department

We repair everything we sell!

1-800-FIX-ENDO

(800-349-3636)

With 38 years of industry experience, ESS stands as the longest-standing and most reliable independent scope service provider, catering to customers globally. How can we assist you?



Save money, and keep your endoscopes at peak performance!

**ScopePlus
Service Contract**

endoscopy.com/Service

Some of the repairs covered:

- Adjust/Replace Control Cables
- Unblock Air/Water Nozzle
- Replace Bending Rubber
- Replace Biopsy Channel
- Replace Light Guide

Save up to **50%** on the cost of a major repair!